


*Your One-Stop-Shop for aseptic Fill
and Finish of biopharmaceuticals*



BIOVIAN

*Contract Development and
Manufacturing of Biopharmaceuticals*

BIOVIAN provides aseptic Fill and Finish, which is the final, critical step in the GMP manufacturing of biopharmaceutical products. Aseptic Fill and Finish is a core competence of Biovian, with an expertise with continuous operation for over a decade. Aseptically filled GMP-certified Drug Products can be used for clinical trials or for commercial purposes.



"When partnering with us there is no need to use several service providers, which can be both time-consuming and expensive."

-Knut Ringbom, CEO

FILL AND FINISH EXPERTISE WITH A BROAD RANGE OF BIOLOGICS

Viral Vectors: AAV, Adenovirus and other BSL2 class Viral Vectors
Plasmid DNA
Recombinant proteins and protein complexes
Dilution buffers
Placebos

FILL AND FINISH SERVICES

Our complete Fill and Finish service chain for Drug Products includes formulation, aseptic Fill and Finish, terminal sterilization or aseptic processing, visual inspection, in-house sterility testing, labeling, packaging and Qualified Person certified release.

ONE-STOP-SHOP SERVICES FOR FILL AND FINISH



FORMULATION



ASEPTIC FILL AND FINISH



VISUAL INSPECTION



LABELING AND PACKAGING




QP CERTIFICATION

ASEPTICALLY FILLED GMP-CERTIFIED DRUG PRODUCT
READY FOR CLINICAL TRIALS OR FOR COMMERCIAL USE

FILL AND FINISH FACILITIES AND STORAGE

Biovian's aseptic filling areas fulfil the requirements for handling hazardous products/biologics in the GMO II group. We also provide intermediate warehousing of the Drug Products at various temperatures (RT, +4°C, -20°C and -80°C) on the same European Medicines Agency (EMA) certified site.



FILL AND FINISH RESOURCES

Dedicated filling lines for BSL1 and BSL2 products
Restricted-Access Barrier System (RABS) equipped filling line for recombinant proteins and plasmid DNA
Range of volumes: from 0,5 mL to 100 mL
Vial-sizes 2R and 10R validated (ask for customized options)
Batch size ranges from 50 to 10 000 vials
BSL1 up to 10 000 vials
BSL2 up to 1000 vials
In-house release testing including sterility testing and 100% visual inspection

*"Our complete Fill and Finish service chain
for Drug Products is part of our core competence."*



*From gene to finished vial –
from the laboratory bench to the clinic*

EXPLORE OUR FILL AND FINISH SERVICES



FOLLOW US ON



CONTACT



contact@biovian.com



Tykistökatu 6 A, Biocity
FI-20520 Turku, Finland



www.biovian.com